


Signal Transduction and Cancer

Generation of High-Quality Data


PD Dr. Ursula Klingmüller
Theodor-Boveri-Group „System Biology of Signal Transduction“

Dynamic Growth and Differentiation Processes




- Qualitative data
- Static description

Kubitzky et al., 2001
Current Biology
Klingmüller et al., 1996
PNAS
Klingmüller et al., 1997
PNAS
Klingmüller et al., 1995
Cell




Dynamic Behavior - Decision Making Process


- Timing
- Amplitude
- Duration
- Spatial Organization




JAK-STAT Signaling Cascade



Mathematical Model (1): JAK-STAT Feed-Forward Cascade




Mathematical Model (2): Cycling STAT-Signaling Pathway




JAK2-STAT5 Pathway

Estimating the dynamical parameters (k)




JAK2-STAT5 Pathway

Quantitative, time-resolved measurements




JAK2-STAT5 Pathway

Fit of Experimental Data: Feed-Forward Cascade (1) versus Cycling (2)




JAK2-STAT5 Pathway


Predicting Steps Most Sensitive for Perturbation



JAK2-STAT5 Pathway

Experimental validation of prediction






Development/ Improvements of Quantitative Immunoblotting


- CCD camera based analysis
- Conversion of arbitrary units to molecules per cell
- Determining steps prone to errors
- Normalizers
- Calibrators

Determination of molecules per cell




- linear detection
- generation of absolute numbers

Determination of the Error of Immunoblotting Measurements




Determining Steps Prone to Error





- uncorrelated error (pipetting, antibody detection) is small
- strong impact of correlated error (gel/transfer inhomogeneities)


Determining Steps Prone to Error



- randomized gel loading destroys correlated error (gel/transfer inhomogeneities)



Time-course of STAT3 Activation in Primary Hepatocytes
Total Cell Lysate




Calibrators

- Prepurification of proteins by immunoprecipitation
- Quantitative analysis of rare proteins from cell lysate by immunoprecipitation
- Error prone multi-step process

Calibrator
(recombinant protein that harbors the epitope recognized by the antibody used for immunoprecipitation)

- added to cell lysates
- used for normalization

Time-course of STAT3 Activation in Primary Hepatocytes
Immunoprecipitation



Generation of High Quality Quantitative Data
Identification of error prone steps - Improvements

Major source for errors:

- SDS-PAGE and transfer

Less important:

- Incubation with antibody and luminol

Improvements:

- Randomized gel loading
- Prepurification by immunoprecipitation
- Use of appropriate normalizers and calibrators

Summary

- New biological knowledge can be generated by data-based mathematical modeling
- High quality data can be generated by quantitative immunoblotting applying randomized gel loading and appropriate normalizers and calibrators

Theodor-Boveri-Group
„Systems Biology of Signal Transduction“

Ursula Klingmüller
Ute Baumann
Achim Heinrich
Andrea Pfeifer
Marcel Schilling
Verena Becker
Peter Nickel
Anna Geist
Sebastian Bohl
Christian Schröter



Olivier Sandra
Ira Swameye

Freiburger Zentrum
für Datenanalyse und
Modellbildung

Jens Timmer
Thorsten Müller
Thomas Maiwald
Markus Kollmann
Kilian Bartholome
Clemens Kreutz